



SONS TRAUMA KIT



MADE IN USA



US PATENTS
PENDING



DEPLOYMENT TRAUMA KIT

PART# DTK

Where do you carry your Trauma Kit...
and what condition is it in?

There's nothing more important than gear readiness when saving a life. Proven in Afghanistan and Iraq, the **Sons Trauma Kit** is engineered for maximum durability and efficiency during trauma treatment. The fit is undetectable, comfortable to carry and instantly accessible.

The contents are packed around injury patterns in order of priority for care and follow all TCCC protocols.

ACCESSIBILITY:

Can be retrieved instantly from under your ballistic protection while in a seated position, in tight quarters, or in confined vehicles. Replaces bulky kits worn are either on the leg, side or on the back. Should the patient need to be dragged out of harms way, this kit will remain on the person. Current kits get torn off in that process.

TECH SPECS:

Weight:	1/2 lbs
Thickness:	1/2 - 3/4 in
Width:	9 in
Length:	12 in
Material:	Medical Grade Vinyl

EASE OF USE:


Designed for quick deployment, the contents are packed around injury patterns in order of priority for care. A full length Ripcord opens the pressure sealed kit.

COMFORT:

Low profile design, lightweight and no nylon pouch needed. Eliminates the danger of getting snagged when exiting tactical vehicles. Does not impede tactical movements, such as crawling, moving through small areas. Offers Buoyancy.

BATTLE TESTED DURABILITY:

Securely fits under a MIL-SPEC Medium, Level III ballistic plate and is protected from blast type attacks. Machine sealed medical grade vinyl and precision RF Welding produces a smooth continuous seam that creates waterproof protection, designed to last 24-30 months, keeping gear protected during long deployments. Outperforms the competition by a ratio of 2 or 3 to 1.




*contents may vary



2" Safety Pins x2



Nitrile Gloves



3" Elastic Bandage



14G Needle



Occlusive Dressing x2*



4" Z-Fold Gauze Bandage x2



Hemostatic Clotting Gauze*



NPA w/ 550*



Casualty Card

RESCUE **E** ESSENTIALS

(866) 711-4843
www.Rescue-Essentials.com