

# RESCUE ESSENTIALS

TACTICAL EMERGENCY MEDICAL SUPPLIES

## MCI BLEEDING CONTROL PACK

The MCI Bleeding Control Pack is a comprehensive mass casualty response bag for active shooter incidents. It provides the ideal Rescue Task Force platform for casualty care in the warm zone of an active incident. Capable of treating up to 8 casualties, the MCI pack is designed for rapid deployment of medical supplies to the scene for treating casualties in multiple locations prior to evacuation to EMS casualty collection. Tourniquets, trauma shears, and gloves are externally mounted for rapid access. Inside the pack are five rip-away kits for rapid deployment to patients.



### Throw & Go Trauma IFAK's (2 Each)

- 2 Pair Nitrile Gloves
- 1 CAT Tourniquet
- 2 Compressed Gauze
- 1 Pressure Dressing
- 1 Hyfin Vented Chest Seal, Twin Pack
- 1 Hypothermia Blanket
- 1 Nasopharyngeal Airway, 28 Fr
- 2 TCCC Casualty Cards
- 1 Sharpie Marker

### External Supplies

- 4 CAT Tourniquets
- 2 Trauma Shears, 7.5"
- 4 SWAT-T Tourniquets
- 1 Glove Dispenser w/ 6 Pairs Nitrile Gloves
- 2 TacMed Blast Bandage
- 1 Rescue Essentials QuikLitter™
- 1 Nylon Webbing, 1"x22.5ft
- 6 TCCC Casualty Cards
- 2 Sharpie Markers

### Bleeding Control Kit

- 6 Compressed Gauze
- 3 Pressure Dressing
- 2 Trauma Dressing, 8" x 10"

### Airway/Chest Kit

- 6 Hyfin Vented Chest Seal, Twin Pack
- 6 Nasopharyngeal Airway, 28 Fr
- 4 ARS Decompression Needle
- 2 Combat Medic Tape
- 4 Pair Nitrile Gloves
- 6 Hand Sanitizers

### Hypothermia Kit

- 6 Hypothermia Blanket



910.830.0286

info@rescue-essentials.com  
8770 Trade St., Leland, NC 28451